

Clinispin™

Micro III

High Speed Micro Veterinary Centrifuge



HIGH SPEED
FOR SMALL VOLUMES

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.



High Speed Micro Veterinary Centrifuge

Clinispin Micro III is a high speed micro veterinary centrifuge equipped with a maintenance free drive, digital display and simple interface for silent and efficient operation. Designed for centrifuging small volumes of biological fluids at speeds of up to 16,500 RPM.

FEATURES

- Rotor allows for quick and easy sample loading
- Maintenance free drive
- Digital display
- Easy to use
- Variable speed setting in increments of 100 RPM up to 16,500 RPM
- Variable timer setting from 1 to 99 minutes
- Displays RPM or RCF (g force)
- Lid lock and imbalance detection safety features



SPECIFICATIONS

Motor Type: Brushless DC Motor
Maximum Speed: 16,500 RPM
Maximum RCF: 18,260xg
Rotor Volume: 12 x 1.5/2.2ml tubes
Timer: 1-99 minutes
Power: AC220-230V 50/60Hz

Weight: 7kg
Dimensions: 300 (L) x 250 (W) x 200 (H) mm
Permitted Temperature Range: 5°C - 40°C
Relative Humidity: ≤80%
Altitude: ≤2000m
Noise (at Maximum RPM): ≤65dB(A)



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

**WOODLEY
EQUIPMENT COMPANY LTD.**