



# Small Multipurpose Veterinary Centrifuge



**VERSATILE.**  
**COMBINATION ROTOR.**

**"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"**

**WOODLEY**  
EQUIPMENT COMPANY LTD.

The Clinispin MPC Combi is a high speed centrifuge equipped with a maintenance free drive, digital display and simple interface for silent and efficient operation. The combination Rotor provides efficient separation of blood and urine samples and accurate Haematocrit determination.

## FEATURES

- 3 pre-defined presets for easy and rapid operation at one press of a button:

- Serum / Plasma (12,000rpm for 2 minutes)
- HCT / QBC® (12,000rpm for 5 minutes)
- Urine (1,000rpm for 1 minute)

Pre-sets can be varied for each centrifugation if required

- Rotor allows for quick and easy sample loading
- Maintenance free drive
- Digital display
- Easy to use
- Quiet operation
- Displays RPM or RCF (g force)



**Includes MPC Combi Rotor for  
18 x HCT/QBC tubes & 4 x 2ml tubes**

## SPECIFICATIONS

- **Motor:** Brushless DC Motor
- **Maximum Speed:** 12,000 RPM
- **Maximum RCF:** 14,119 x g (HCT / QBC® tubes)  
9,982 x g (2 ml tubes)
- **Maximum Time:** 15 minutes
- **Weight:** 6kg
- **Dimensions:** 365 x 291 x 134 mm (L x W x H)
- **Maximum Volume:** 18 x 75mm (HCT/QBC®)  
4 x 2ml (Micro Tube)
- **Electrical Requirements:** 24V, 6.25 A



+44 (0) 1204 669033 (Option 1)

sales@woodleyequipment.com

www.woodleyequipment.com

Old Station Park Buildings  
St. John Street  
Horwich  
Bolton  
Lancashire  
BL6 7NY  
UK

**"PASSIONATE ABOUT  
LABORATORY DIAGNOSTICS"**

**WOODLEY**  
EQUIPMENT COMPANY LTD.