

# Clinitek Status®+ Urine Analyser



Available for hire or purchase for your clinical trial.

## Overview

The Clinitek Status®+ analyser is easy to use and offers the broadest menu available on any Point Of Care automated urine analyser. In addition, this urine analyser features new automatic checks (Auto-Checks\*).

The analyser automatically checks each test strip for humidity exposure, common sample interferences and strip identification for Siemens test strips. Together, these provide improved clinical information.

## Operating Feature

- Easy to Use and Convenient
- Simple, intuitive touch screen operation
- Load test – the analyser does the rest, automatically times and reads every specimen
- Fast results – in-office test results in about 1 minute
- Automatic printed report – no manual transcription needed Improved Data Integrity
- Analyser along with patented test strip technology enables Auto-Checks
- Detects common sample interferences and makes a printed note of any potentially compromised result\*\*
- Checks each strip for humidity to prevent reporting of falsely elevated results
- Auto-identification of test strips, no need for operators to key in manually
- Removes subjectivity of visually read results

\*Only available when using Siemens IR or colour bands

\*\*Feature not available in the U.S.

**Please contact us for a quotation for the hire or purchase of this product**

### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: [hello@woodleytrialsolutions.com](mailto:hello@woodleytrialsolutions.com)

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

## Specification

- Routine Urine Tests: albumin, bilirubin, creatinine, glucose, ketone, leukocytes, nitrite, pH, protein, specific gravity and urobilinogen
- Kidney Disease Test: albumin-to-creatinine ratio
- Pregnancy Test: human Chorionic Gonadotropin (hCG)
- Computer Interface: uni-directional via serial port (RS232)
- Instrument Memory: 950 patient test results, 700 operator IDs
- Interface: touch screen display
- Power Requirements: 100 - 240V AC, 50-60Hz (with in-line lead), or 6 x AA batteries
- Operating Temperature: 18°C - 30°C
- Operating Humidity: 18 - 80% RH
- Calibration: automatic, self-calibrating
- Weight: 1.6 kg
- Dimensions: 158h x 171w x 272d mm



Product description : Clinitek Status®+ Urine Analyser

Product specifications can be subject to change without prior notice. Version: V2 10/19

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company