

Clinitest® hCG Pregnancy Test

Available for hire or purchase for your clinical trial.



Overview

The Clinitest® hCG pregnancy test allows for simple and reliable testing using the Clinitek Status®+ analyser.

The Clinitest® hCG pregnancy test gives you fast, easy instrument-read-and-reported test results - taking the question out of results and the interpretation out of your hands. Just apply the patient's sample and let the Clinitek Status® analyser do the rest.

Operating Feature

- Quick, reliable instrument-read results eliminate interpretation and provide accurate results for your patient
- The Clinitek Status Analyser precisely times the reading of results, eliminating the hands-on waiting time
- Printed results eliminate transcription errors and save time
- CLIA-waived
- Reduced user involvement and time
- Enhanced uniform test method
- Decreased user uncertainty
- Less opportunity for error
- More objective results
- Increased time for patient care
- Consistent, reliable results every time

Specification

- Quality Control: two procedural controls
- Storage: room temperature
- Sample Size: 200µL urine for a single test
- Time To Result: 5 minutes or less to confirm result
- Method: lateral flow, chromatographic result
- Clinical Accuracy: >99%
- Sensitivity: equal or greater than 25 ml U/ml



Please contact us for a quotation for the hire or purchase of this product

Product description - Clinitest® hCG Pregnancy Test

Product code

Product specifications can be subject to change without prior notice.

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: enquiries@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company