

CoaguChek® XS Plus System

Available for hire or purchase for your clinical trial.



Overview

The CoaguChek® XS Plus coagulation monitoring system will help you to provide on the spot high quality PT/INR results, allowing you and your patients more quality time and increased peace of mind.

The CoaguChek® XS Plus has several innovations; easy portability, top and side blood dosing options, extensive data management options and safety features, e.g. onboard quality control (QC) for strip integrity check, optional QC lockout to meet potential regulatory requirements and more.

Operating Feature

- Full graphical touch screen for simple and fast handling
- Audit facility with patient and user identification
- Download capabilities and integrated data management functions
- Accurate results within 1 minute
- 500 result memory
- Mains and battery power supply for in-clinic or home visit use
- CoaguChek® XS Plus ProthrombinTime (PT) liquid controls available
- Detection System: amperometric (electrochemical) determination of the PT/INR after activation of the coagulation with human recombinant thromboplastin
- QC lockout, patient & user ID, protection with administrator ID

Specification

- Sample Type: fresh capillary whole blood or non anti-coagulated venous blood
- Sample Size: 10 µl
- Operating Temperature: 18°C - 32°C
- Humidity and Height: 10 - 85% RH - 4,300 m
- Measuring Range: INR: 0.8 - 8.0
- Memory: 500 test results with date and time, 60 code chip records
- Power Options: 100 - 240V, 50 - 60 Hz, 400mA or 4 x 1.5V AA batteries
- Battery Life: approx. 80 tests
- Weight: 311g without batteries
- Dimensions: 185h x 98w x 42d mm



Please contact us for a quotation for the hire or purchase of this product

Product description : CoaguChek® XS Plus System
Product specifications can be subject to change without prior notice. Version: V2 10/19

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company