

Counter Height Medical Laboratory Pharmacy Freezer 5.0 cu.ft./141 L

Available for hire or purchase for your clinical trial.



Overview

This counter height medical laboratory freezer has excellent temperature stability and management and is suitable for use in medical, scientific, pharmaceutical and laboratory facilities where precise temperature conditions are required.

The 5.0 cu. ft. 22" wide counter height freezer is designed for free standing use. Two open shelves include cooling coils for optimum performance, with additional door racks for easy storage.

Operating Features

- Slim counter height dimensions
- Manual defrost
- Digital thermostat
- High temperature alarm with an external temperature readout on the door
- Hospital grade cord with green dot plug
- Door storage
- Side mounted keyed lock
- Two removable wire baskets
- Fully finished white cabinet
- Flat door liner
- Access port - unit is supplied with factory installed port for users to connect monitoring equipment
- 2 Level legs
- US Electrical Safety - UL
- Canadian Electrical Safety - ULC
- 100% CFC Free

Specification

- Capacity: 5.0 Cu. Ft. / 141 Litres
- Temperature Range: - 18 to -12 °C / 0 to 10°F
- Amps: 1.1
- Input Voltage: 115 V AC/60 Hz
- Weight: 90 lbs
- Dimensions External: 33.13"h x 23.0"w x 25.75"d
- Dimensions Internal: 29.5" h x 18.63"w x 19.0"d



Please contact us for a quotation for the hire or purchase of this product

Product description: Counter Height Medical Laboratory Freezer 5.0 Cu. Ft.

Product Code: WD7FRZ089

Version: V2 10/19

Product specifications can be subject to

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company