

# Counter Top Medical Laboratory Freezer 1.4 Cu.Ft./40 L

Available for hire or purchase for your clinical trial.



## Overview

This compact counter top medical laboratory freezer has a stabilized temperature performance to be used in medical, scientific, pharmaceutical and laboratory facilities.

The 1.4 cu. ft. manual defrost interior offers static cooling for stable results and a removable wire shelf allows easy storage options.

With high and low temperature audible alarm that sounds if the interior temperature rises or falls out of set range. The freezer has external temperature readout.

## Operating Features

- Compact dimensions ideal for counter top use
- Manual defrost
- Digital thermostat
- Factory installed lock
- High/low temperature alarm
- NIST calibrated temperature display
- High grade cord with green dot plug
- One removable wire basket
- One piece interior liner
- Flat door liner
- Access port - unit is supplied with factory installed port for users to connect monitoring equipment
- 2 Level legs
- US Electrical Safety - ETL
- Canadian Electrical Safety - ETL-C
- 100% CFC Free

## Specification

- Capacity: 1.4 Cu. Ft. / 40 Litres
- Temperature: - 20 to -15 °C / - 5 to +5 °F
- Amps: 0.8
- Input Voltage: 115 V AC/60 Hz
- Weight: 35 lbs
- Dimensions External:  
21.38 x 18.5 x 17.63 inches H x W x D
- Dimensions Internal:  
14.38 x 13.5 x 13.25 inches H x W x D



**Please contact us for a quotation for the hire or purchase of this product**

Product description: - Counter Top Medical Laboratory Freezer 1.4 Cu. Ft. / 40 Litre  
Product Code: WD7FRZ088  
Version: V2 10/19  
Product specifications can be subject to

### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company