

# CT321 PC based ECG Machine

Available for hire or purchase for your clinical trial.



## Overview

The CT321 is a PC based interpretive ECG Machine with Bluetooth connectivity. The innovative Medical Bluetooth Technology removes the need for a USB cable connection, allowing greater flexibility in large clinic rooms.

The CT321 acquires the ECG on any computer/laptop within your healthcare network, allowing you to directly transfer the patient's ECG report into your Electronic Medical Records (EMR). It is compatible with all leading EMR systems.

## Operating Features

- Supplied complete with patient cable, alligator clips, 500 disposable electrodes, cardio software (no licence fees), Bluetooth dongle, 2 rechargeable lithiumion batteries and battery charger
- Cardio software automatically exports a pdf ECG file to a selected destination
- Easy attachment of the EMR to email to specialist ECG centres and for use in tele-health applications
- Cardio ECG software functions: single beat analysis, cumulative complexes and measurement table
- Paperless ECG records, if hard copy is required, the ECG can be printed in a choice of formats complete with grid onto plain A4 copier paper
- Single patient database viewed throughout practice network
- Advanced software filters for baseline, electrical interference, anti drift and muscle tremor interference, providing clear traces
- Accessories: Snapclip Connectors (Pack x 10)
  - Q Trace Electrodes
  - Alligator clips 4mm (Pack x 10)

## Specification

- Digital pacemaker detection
- Sampling rate of up to 4000 Hz
- Robust defibrillator proof unit
- Impedance measurement: all electrodes highlighted with automatic quality indication to ensure the perfect ECG
- Size: 16 x 8.5 x 2.5cm
- Weight: 300g



**Please contact us for a quotation for the hire or purchase of this product**

Product description: CT321 PC Based ECG Machine

Product Code: WD7ECG013

Version: V2 10/19

Product specifications can be subject to change without prior notice

### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company