

# CT6i Interpretive ECG Machine

Available for hire or purchase for your clinical trial.



## Overview

The CT6i ECG machine features user friendly controls and settings, and comes complete with an interpretation and measurements program.

The A4 presentation means you can print either all 12 leads per page or 4 x 3 leads plus rhythm strip, intervals and interpretation on one A4 page for easy photocopying and filing of records.

Interpretation and measurements included.

## Operating Features

- 12 channel, A4 presentation, interpretive ECG
- Copy button for extra traces
- Useful "Leads Off" warning with lead identification
- A4 print out with a choice of up to 12 leads per page or 4 x 3 leads plus rhythm strip, intervals and interpretation
- Auto-start one touch automatic operation or manual control if preferred
- 12 lead simultaneous acquisition
- Mains, base line and muscle filters, for clear traces
- Interpretation report, with or without complexes, measurements and rhythm strips can be preset for your preferred one touch operation

## Specification

- Measurement Table: RR, P, PR, QRS, QT and QTC intervals, and the P, QRS and T axis
- Mains and battery operation with automatic recharge
- Fully charged battery will provide up to 400 ECG's
- Weight: 5kg with battery installed
- Dimensions: 104h x 399w x 338d mm
- Supplied complete with mains cable, patient cable, pack of A4 paper, 500 disposable electrodes and alligator clips.



**Please contact us for a quotation for the hire or purchase of this product**

Product description: CT6i ECG Machine

Product code - WD7ECG002

Version: V2 10/19



Product specifications can be subject to change without prior notice

### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: [hello@woodleytrialsolutions.com](mailto:hello@woodleytrialsolutions.com)

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company