

CT8000P Interpretive ECG Machine

Available for hire or purchase for your clinical trial.



Overview

The CT8000P ECG machine features user friendly controls and settings and comes complete with an Interpretation and Measurements program. It features an LCD screen and alpha-numeric keyboard.

Equipped to provide paperless ECG options when combined with Achimed Passport (Open) Software.

Operating Features

- 12 channel high quality interpretive ECG with LCD screen and alpha-numeric keyboard
- High quality adult and paediatric interpretive software
- Easy selection of traces from 12 x 1, 6 x 2 or the popular 4 x 3 plus Rhythm Strip one page report
- On-screen lead test and lead-off indication
- One touch automatic or instant real time operation
- High quality digital filters for clear traces (mains, base line and muscle tremor)
- Internal memory for up to 45 stored ECG records
- All parameters including filters are easily programmed through the on-screen menu for a one touch ECG
- Dedicated one touch keys for storing, printing and transmitting
- Easy-load paper tray (220 A4 pages)
- Interpretation report, with or without complexes, measurements and Rhythm Strips can be preset for your preferred one touch operation
- Mains and battery operation

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office


t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

Specification

- Measurement Table: RR, P, PR, QRS, QT and QTC intervals, and the P, QRS and T axis
- Mains and battery operation
- Print Speed: 5, 10, 25, 50 mm/sec
- Weight: 8 lbs/ 5kg (with paper)
- Dimensions: 13h x 15w x 13d inches / 330h x 390w x 330d mm
- Supplied complete with mains cable, patient cable, pack of A4 paper, 500 disposable electrodes and alligator clips.
- Upgradeable to Spirometry, Stress Test and Thrombolysis Indication, please enquire for full details. Spirometry includes pre- and postmedication with results displayed and printed with diagnosis statements based on international standards



Please contact us for a quotation for the hire or purchase of this product

Product description: CT8000P ECG Machine

Product code - WD7ECG003

Version: V3 10/19



Product specifications can be subject to change without prior notice

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company