

CT320 PC based ECG Machine

Available for hire or purchase for your clinical trial.



Overview

The CT320 is a PC based interpretive ECG Machine with USB connectivity. The CT320 acquires the ECG on any computer/laptop within your healthcare network, allowing you to directly transfer the patient's ECG report into your Electronic Medical Records (EMR). It is compatible with all leading EMR systems.

Cardio software automatically exports a pdf ECG file to a destination of your choice enabling ease of attachment to the EMR or for emailing to specialist ECG centres and use in tele-health applications.

Operating Features

- Easy to use interface
- Supplied complete with USB connection cable, full length patient cable, alligator clips, 500 disposable electrodes and cardio software (no licence fee)
- Cardio software automatically exports a pdf ECG file to a selected destination
- Cardio ECG software functions: single beat analysis, cumulative complexes and measurement table
- Paperless ECG records, if hard copy is required, the ECG can be printed in a choice of formats complete with grid onto plain A4 copier paper
- Single patient database viewed throughout practice network
- Accessories:
 - Snapclip Connectors (Pack x 10)
 - Q Trace Electrodes
 - Alligator clips 4mm (Pack x 10)

Specification

- Digital pacemaker detection
- Sampling rate of up to 4000 Hz
- Robust defibrillator proof unit
- Power supply via USB connection cable
- Size: 16 x 8.5 x 2.5cm
- Weight: 300g



Please contact us for a quotation for the hire or purchase of this product

Product description: CT320 PC Based ECG Machine

Product Code: WD7ECG014

Version: V2 10/19

Product specifications can be subject to change without prior notice



UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company