

Digital Class III Flat Scales

Available for hire or purchase for your clinical trial.



Overview

It is the Department of Health recommendation that all new scales purchased for the purpose of weighing patients should as a minimum achieve a class III standard.

These scales have been developed with this in mind. They are a low cost alternative scale which comply with the regulations but still achieve on performance.

These compact personal scales with their durability, functionality and size permit medically precise weighing results in even very confined spaces.

Operating Features

- Automatic power down function
- Lightweight
- Outstanding surface hardness and impact strength
- 4 leveling bases
- Slip proof surface
- Easy to clean
- Large levelling base
- Slip proof standing surface
- Tip-on automatic function

Specification

- Capacity: 200 kg
- Scale division: 100g < 150 kg > 200g
- Graduation: 200 g
- Power supply: batteries
- Approval: CE approval class (III)
- Weight: 3.6 kgs
- Dimensions: 59h x 321w x 356d mm



Please contact us for a quotation for the hire or purchase of this product

Product description: - - Digital Class III Flat Scales

Product code - WD7PAT006

Product specifications can be subject to change without prior notice. Version: V2 10/19

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company