

Digital Column Scales

Available for hire or purchase for your clinical trial.



Overview

These column scales are versatile, low maintenance and site-independent with economic battery operation. They are ideal for daily use in hospitals and doctors practices.

They are equipped with several intelligent functions - for example the BMI function permits a reliable evaluation of the nutritional condition of the patient.

They feature transport castors so they are easy to move around.

Operating Features

- Functions:
 - BMI
 - HOLD function
 - TARE function
 - lbs / kg switch over
 - automatic / manual switch off
- Smooth running transport castors
- Low maintenance
- Site-independent
- Battery operated
- Energy-saving and cost-effective
- Main adapter (optional)

Specification

- Capacity: 450 lbs
- Graduation: 0.2lbs
- Graduation measuring range: 1/8 inch
- Power supply: batteries / main adapter (optional)
- Measuring rod (supplied)
- Measuring range 24 to 78 inch
- Weight: 13.7 lbs
- Platform Dimensions inches:
H 3.0" x W 10.7" x D 11"
- Dimensions inches:
H 32.7" x W 11.6" x D 16.4"



Please contact us for a quotation for the hire or purchase of this product

Product description: - Digital Column Scales

Product code - WD7PAT092

Product specifications can be subject to change without prior notice. Version: V3 10/19

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company