

► Digital Column Scales

available for hire or purchase for your clinical trial.

Overview

These column scales are versatile, low maintenance and site-independent with economic battery operation. They are ideal for daily use in hospitals and doctors practices.

They feature transport castors so they are easy to move around.

They are equipped with several intelligent functions - for example the BMI function permits a reliable evaluation of the nutritional condition of the patient.

Specification

- Capacity: 200 kg
- Graduation: below 150 kg - 100g increments
above 150 kg - 200g increments
- Power supply: batteries / main adapter (optional)
- CE approval class: (III)
- Measuring rod (optional)
- Weight: 6.3 kg
- Platform Dimensions: 75h x 272w x 280d mm
- Dimensions: 831h x 294w x 417d mm

Operating Features

- Functions:
 - BMI
 - HOLD function
 - TARE function
 - lbs / kg switch over
 - automatic / manual switch off
- Smooth running transport castors
- Low maintenance
- Site-independent
- Battery operated
- Energy-saving and cost-effective
- Main adapter (optional)



Please contact us for a quotation for the hire or purchase of this product.

Product description - Digital Column Scales
Product code - WD7PAT004



UK Office

T: +44(0)8456 777001

F: +44(0)8456 777002

E: enquiries@woodleyequipment.com

USA Office

T: 1 800 471 9200 / 1 718 606 0516

F: 1 718 606 0955

www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.