

Digital Flat Scales

Available for hire or purchase for your clinical trial.



Overview

A 6.5 foot cable connects the display with the weighing base. Advantage: the display can be set up separately from the base and read from its position on the wall or desk. Additional convenience is supplied by the many extra weighing functions: BMI for nutritional conditions, HOLD to keep measurements displayed and TARE for weighing small children held in an adult's arms. It is suitable for stationary or mobile use thanks to its handle, compact size, low weight and a leveling base for set-up on all types of surfaces.

Operating Features

- Automatic switch-off
- Stylish and sleek design
- BMI function providing nutritional information
- HOLD
- TARE
- 6 x Type AA 1.5V batteries
- Mobile use

Specification

- Capacity: 550 lbs
- Graduation: 0.2lbs < 330lbs > 0.5lbs
- Power Supply: batteries
- Weight: 9.9lbs
- Dimensions: product: inches: H 2.4" x W 12.6" x D 14"
- Dimensions: platform inches: H 2.4" x W 11.3" x D 10.9"
- Cable Length: 6.5 feet
- Functions: HOLD, BMI, Automatic switch-off, Mobile measuring, TARE



Please contact us for a quotation for the hire or purchase of this product

Product description: Digital Flat Scales (869)

Product code - WD7PAT093

Product specifications can be subject to change without prior notice. Version: V2 10/19

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company