

Flat Scales For Mobile Use

Available for hire or purchase for your clinical trial.



Overview

These scales are suitable for mobile use. Lightweight, stable, compact and easy to clean they are equipped with a Mother-Child function.

The slip proof standing surface and steel element in the base allows for patients up to 200 kg to be weighed easily. The scales are easy to use with the large display and operating element integrated into the base.

Operating Features

- Mother - Child function for trouble-free weighing of small children
- Functions: Mother/child function, tip-on automatic switch on, automatic switch off
- Large levelling base for high stability
- Easy to transport - scale weighs 4.2 kg
- Flat construction and slip proof standing surface
- Battery powered for use in locations without electrical supply

Specification

- Capacity: 200 kg
- Graduation: 100 g < 150 Kg > 200 g
- Power Supply: batteries
- CE Approval Cass: (III)
- Weight: 4.2 kg
- Dimensions: display: 60h x 321w x 356d mm



Please contact us for a quotation for the hire or purchase of this product

Product description: - Flat Scales for Mobile Use (877

Product code - WD7PAT007

Product specifications can be subject to change without prior notice. Version: V2 10/19

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company