

Drugs Of Abuse Test Panels

Available for hire or purchase for your clinical trial.



Overview

The Triage® drugs of abuse panel is a small, handheld urine test that provides qualitative results for up to 9 commonly abused drugs in about 10 minutes. The range of panels provide an extensive range of testing.

Amphetamine (AMP), Methamphetamine (mAMP), Barbiturates (BAR), Benzodiazepines (BZO), Cannabinoids (THC), Cocaine (COC), Opiates (OPI), PCP (PCP), Tricyclic antidepressants (TCA), Methadone (MTD), Propoxyphene (PPX) and Acetaminophen.

Operating Feature

- Excellent discrimination around assay cut offs
- Ensures accurate patient and proficiency results
- Provide superior specificity
- Minimise false negative results
- Template permits hard copy archiving of results
- Greater cost efficiency
- Accepted by COLA, JCAHO and CAP
- Positive specimens produce a clear visual result
- Negative specimens do not produce a visual result

Please contact us for a quotation for the hire or purchase of this product

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

Specification

- Tests a wide range of analytes -

TOX Drug Screen

AMP, mAMP, BAR, BZO, THC, COC, OPI, PCP, TCA, MTD, Acetaminophene

DOA Panel

AMP, mAMP, BAR, BZO, THC, COC, OPI, PCP

DOA Panel + TCA

AMP, mAMP, BAR, BZO, THC, COC, OPI, PCP, TCA

DOA Panel + MTD

AMP, mAMP, BAR, BZO, THC, COC, OPI, TCA, MTD

DOA Panel + PPX

AMP, mAMP, BAR, BZO, THC, COC, OPI, PCP, TCA, PPX

Express Test 5

AMP, mAMP, THC, COC, OPI, PCP



Drugs Of Abuse Test Panels

Product specifications can be subject to change without prior notice. Version: V2 10/19

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company