

www.woodleyequipment.com

WOODLEY EQUIPMENT / QUANTUM

in association with

RANDOX LABORATORIES LIMITED

introduce

ECLIPSE

CLINICAL CHEMISTRY ANALYSER



Key Benefits

- 90 tests/hour
- Small sample requirement - minimum volume 120ul
- Convenient onboard refrigerated reagent storage
- Flexible test combinations - full profile or single analyte
- Auto re run when results outside technical range of system
- Up to 50 species specific reference ranges can be set
- Minimal user maintenance

Tel: +44 (0)1204 669033

Fax: +44(0)1204 669034

Email: sales@woodleyequipment.com

Web: www.woodleyequipment.com

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK

WOODLEY
EQUIPMENT COMPANY LTD.

Now Introducing
The
ECLIPSE

FULLY AUTOMATED VET CHEMISTRY ANALYSER

Randox Laboratories have selected
Woodley Equipment/Quantum Vet Diagnostics to be their
exclusive Veterinary representative in the UK and Ireland.

Uniquely Structural Pricing

Woodley / Quantum Vet Diagnostics are fully aware of the financial constraints that some practices find themselves in and offer an **"All Inclusive"** approach to the pricing structure, which not only includes a regular supply of reagents, regular maintenance visits, all parts and servicing but also, to help assist the practice in their long-term financial planning, the price is fixed for five years.

The Woodley / Quantum System Includes As Standard:

- An All Inclusive Reagent, Support & Maintenance Package
- Out Of Hours Technical Support
- Fixed Price for 5 Years
- Regular Preventative Maintenance Visits
- Part of the Woodley/Quantum Complete Laboratory System including Biochemistry, Haematology, Blood Gases & Coagulation
- Optional Endocrinology Analyser (TOSOH AIA-360)

Customer defined profiles and/or tests can be programmed from the following list:

Albumin, Alkaline Phosphatase, ALT, AST, Amylase, Bilirubin (Total), Calcium, Cholesterol, CK, Creatinine, Glucose, GGT, Phosphate, Protein (Total) and Urea.

Further tests which may also be added to the Woodley/Quantum Series of Practice Analysers, are available, including:

BHB, Bile Acids, Copper, Fibrinogen, Fructosamine, GLDH, GSH-Px (as Selenium), IgG, LDH, Lipase, Magnesium, SAA, T4 (Total Thyroxin), Triglycerides, Uric Acid, Urine Creatinine, Urine Protein, Zinc and Phenobarbitone.

Tel: +44 (0)1204 669033

Fax: +44(0)1204 669034

Email: sales@woodleyequipment.com

Web: www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK