

# Eco Temp Basic Digital Thermometer



Available for hire or purchase for your clinical trial.

## Overview

The Eco Temp Basic Digital Thermometer has a round tip designed for the patient to use safely and conveniently during all measurement methods.

The Eco Temp Basic is water resistant and has an exchangeable battery.

The on/off button is easy to use.

For oral, axillary and rectal use.

## Operating Features

- Digital fever thermometer for oral, axillary and rectal use
- Normal response in 60 seconds
- Water resistant
- 3 digit display (1 decimal place)
- Reading display in Celsius or Fahrenheit
- Beeper signal when measurement is finished
- Memory for last reading
- Automatic switch off saves on battery-life
- Supplied with a protective case
- Exchangeable battery

## Specification

- Measurement Range: 32 to 42°C
- Accuracy of measurement +/- 0.1°C
- Power Supply: 1.5V DC, 1 LR41 Alkaline-Magnesium Button Battery
- Power Consumption: 0.1 mW
- Battery Life: Approximately 2 years or more (3 times a day)
- Temperature Display: 3 digit, +°C display in 0.1 degree increments
- Operating Temperature: +10°C to +40°C
- Relative Humidity: 30 - 85% RH
- Weight: approximately 12g including battery)
- Dimensions: 19.4 x 132.5 x 10 mm



**Please contact us for a quotation for the hire or purchase of this product**

Product description: - Eco Temp Basic Digital Thermometer

Product Code:

Product specifications can be subject to change without prior notice. Version: V2 10/19

### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: [hello@woodleytrialsolutions.com](mailto:hello@woodleytrialsolutions.com)

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company