

EZ-SmartTemp™ Min-Max Thermometer

Available for hire or purchase for your clinical trial.



Overview

The EZ-SmartTemp™ Min-Max thermometer attaches to the exterior of the unit and measures the internal temperature via a probe. The probe is sealed in a miniature bottle filled with non toxic glycol.

The triple display simultaneously shows the high, low and current temperature. Alarm feature alerts if the temperature rises above or falls below a set point.

The EZ-SmartTemp™ Min-Max thermometer is supplied with a calibration certificate.

Operating Features

- Digital display making it easy to read
- Small compact design
- Easy to use
- Displays minimum, maximum and current temperature
- Ideal for a variety of temperature monitoring solutions
- Supplied with a calibration certificate
- The alarm is programmable in 1°C increments
- The probe is sealed in a 5ml bottle. The bottle mimics vaccine temperature and sits in a tray
- 3 mm micro cable permits doors to close on it

Specification

- Temperature Range: -50°C to +70°C
- Accuracy: ± 1°C
- Resolution: ± 1°C / °F
- Battery: 1.5V AAA
- Battery Life: 1000 hours
- Sensor Type: Thermistor
- Display: Custom dual LCD
- Weight: 113g
- Dimensions: 108h x 70w x 20d mm



Please contact us for a quotation for the hire or purchase of this product

Product description: - EZ-SmartTemp™ Min Max Thermometer

Product Code: WD7TEM001

Product specifications can be subject to change without prior notice. Version: V2 10/19

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company