

Flex Temp Smart Digital Thermometer



Available for hire or purchase for your clinical trial.

Overview

The Flex Temp Smart patient thermometer is a high quality, elegant and functional digital thermometer, developed for use in various settings.

This thermometer has a special flexible tip that ensures maximum user comfort during all measurement methods. The compact design is water resistant and the large battery capacity ensures a long life time and operation.

For oral, axillary and rectal use.

Operating Features

- Fast 10 seconds measurement (Rectal: approx. 10 secs, Oral: approx. 20 secs, Armpit: approx. 25 secs.)
- Battery life is approximately 2 years
- Flexible tip for maximum comfort during measurement
- Water resistant
- Memory function recalls last reading
- Temperature can be displayed in Celsius and Fahrenheit
- Measurement method is actual and non predictive
- Beeper guided operation
- Weight is approximately 12g (with battery installed)
- External dimensions 19.4(w) x 132.5 (l) x 10.0(d) mm
- Easy to use placement of the ON/OFF button
- Easy to read display 3 digit (°C), or 4 digit (°F)
- Protective Case included
- 15 language instruction manual

Specification

- Measurement Range: 32.0°C to 42.0°C (89.6°F to 107.6°F)
- Measurement Accuracy: $\pm 0.1^{\circ}\text{C}$ ($\pm 0.2^{\circ}\text{F}$)
- Voltage: 1.5V DC, 1 LR41 Alkaline-Manganese Button Battery
- Power Consumption: 0.1 mW
- Operating Temperature: +10°C to +40°C (+50° to +104°F)
- Operating Humidity: 30 to 85% RH
- Sensing Unit: Thermistor



Please contact us for a quotation for the hire or purchase of this product

Product description: - Flex Temp Smart Digital Thermometer

Product code WD7PAT014

Product specifications can be subject to change without prior notice. Version: V2 10/19

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company