

Gilson Pipette Classic P5000 1-5ml

Available for hire or purchase for your clinical trial.



Overview

The Gilson Pipette Classic P5000 model is robust, accurate and precise pipette. It is ideal for large-volume dispensing of aqueous fluids of moderate viscosity and density. This model is equipped with a polyethylene filter that prolongs piston life and prevents contamination.

Gilson Pipette Classic P5000 model is a fully adjustable, air-displacement pipette with the selected volume shown on the digital indicator (Volumeter).

Please note that this model does not have a tip ejector.

Specification

- | | | | |
|--------------------------|-----------|------|-----------|
| • Volume (µL): | Min. 1000 | 2000 | Max. 5000 |
| • Systematic Error (µL): | ± 12, | ± 12 | ± 30 |
| • Random Error (µL): | ≤ 3 | ≤ 5 | ≤ 8 |



Please contact us for a quotation for the hire or purchase of this product

Product description - Gilson Pipette Classic P5000 1-5ml

Product Code: WD7LAB008

Version: V2 10/19

Product specifications can be subject to change without prior notice.

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company