

# Gilson Pipette Classic P5000 1-5ml

Available for hire or purchase for your clinical trial.



## Overview

The Gilson Pipette Classic P5000 model is robust, accurate and precise pipette. It is ideal for large-volume dispensing of aqueous fluids of moderate viscosity and density. This model is equipped with a polyethylene filter that prolongs piston life and prevents contamination.

Gilson Pipette Classic P5000 model is a fully adjustable, air-displacement pipette with the selected volume shown on the digital indicator (Volumeter).

Please note that this model does not have a tip ejector.

## Specification

• Volume (μL):	Min. 1000	2000	Max. 5000
• Systematic Error (μL):	± 12,	± 12	± 30
• Random Error (μL):	≤ 3	≤ 5	≤ 8



**Please contact us for a quotation for the hire or purchase of this product**

Product description: Gilson Pipette Classic P5000 1-5ml

Product Code: WD7MISC195 (US)

Version: V2 10/19

Product specifications can be subject to change without prior notice

### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: [hello@woodleytrialsolutions.com](mailto:hello@woodleytrialsolutions.com)

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company