

InSight™

Veterinary HCT Meter

The InSight™ HCT meter is a low cost handheld meter for the measurement of canine, feline and equine haematocrit and total haemoglobin in capillary and venous whole blood.

- Calibrated for cats, dogs and horses
- Accurate and precise
- Handheld
- 30 second test time
- One drop whole blood
(Venous, Capillary or EDTA)
- Stores 250 results
- HCT results reported as %
- Room temperature strip storage
- LCD Screen



Calibrated and validated specifically for use in cats, dogs and horses

Tel: +44(0)1204 669033 (Option 1)
Fax: +44(0)1204 669034
Email: sales@woodleyequipment.com
Web: www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK

InSight™

HCT Meter



Kit Contents

- InSight™ HCT Meter
- InSight™ HCT Test Strips (10)
- InSight™ HCT Code Cards – Cat, Dog & Horse
- Lancets (10)
- Lancing Device
- Instruction Manual
- CR2032 Batteries (2)
- Carry Case

Technical Specifications

Parameters:	Haematocrit and Haemoglobin
Measuring Principle:	Optical reflectance
Strip:	InSight™ HCT Test Strips
Display:	LCD
Weight:	<58g (Batteries included)
Outer Dimensions:	102mm x 50mm x 19mm
Power:	DC 6V (Two CR2032 Batteries)
Operating Temperature:	15°C - 30°C

Tel: +44(0)1204 669033 (Option 1)
Fax: +44(0)1204 669034
Email: sales@woodleyequipment.com
Web: www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK