

Height Measure

Available for hire or purchase for your clinical trial.



Overview

This mobile height measure is distinguished by its top-quality materials and a cleverly designed assembly system.

The stability of this height measure is unmatched. It is particularly suited for use in doctors' practices and hospitals.

When disassembled it can be conveniently transported to any location.

Operating Features

- Fast and simple assembly
- Collapsible measuring rod
- Robust base plate guarantees a sturdy stance
- Base has no raised sides or edges
- Wall spacer to prevent wobbling and trembling
- Large and stable head piece for accurate results
- Made from firm non-warping plastic
- Two sided scale on the measuring rod
- Can be dismantled in just a few easy steps
- Adapter element for flat scales (optional)

Specification

- Measuring Range: 20 - 205 cm / 8 - 81"
- Graduation: 1 mm / 1/8"
- Weight: 3.6 kg
- Dimensions: 2145h x 328w x 574d mm



Please contact us for a quotation for the hire or purchase of this product

Product description: -- Height Measure (217)

Product code - WD7PAT010

Product specifications can be subject to change without prior notice. Version: V3 10/19

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company