

# HemoCue® WBC

## Point of Care Testing

For Total White Blood Count



**SIMPLE TO USE.**

**RESULTS IN 2 MINUTES.**

**TOTAL WBC IN BLOOD/SYNOVIAL FLUIDS**



**"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"**

**WOODLEY**  
EQUIPMENT COMPANY LTD.

# HemoCue® WBC

## Point of Care Testing

HemoCue WBC Analyser provides a total White Blood Count.

**Simple to use and results in 2 minutes.**

### KEY FEATURES

- WBC count available using whole blood and synovial fluid
- Small sample volume - 10µl
- Immediate results using proven technology
- Capillary, venous whole blood and synovial fluid can be used
- Factory calibrated
- Powered by either battery or AC adapter for portability

### WHY CHOOSE HEMOCUE?

- **Immediate Results** - test results within 2 minutes
- **Simple to Use** - saves time, work and money
- **Maximum Flexibility** - portable and lightweight
- **Stable Measurement Values** - correct values for correct treatment

**HemoCue WBC  
+  
epoc  
=  
Rapid  
Pathological  
Screening  
24 hours  
a day**



**HemoCue WBC  
+  
InSight HCT  
Meter  
=  
Rapid  
Haematological  
Assessment**



### EASY TO USE IN 3 SIMPLE STEPS



1. The microcuvette is automatically filled with the correct amount of blood/fluid



2. Place the microcuvette in the analyser



3. Read the value on the digital screen within 2 minutes

+44 (0) 1204 669033 (Option 1)

sales@woodleyequipment.com

www.woodleyequipment.com

Old Station Park Buildings  
St. John Street  
Horwich  
Bolton  
Lancashire  
BL6 7NY  
UK

**"PASSIONATE ABOUT  
LABORATORY DIAGNOSTICS"**

**WOODLEY  
EQUIPMENT COMPANY LTD.**