

i-Smart 30 VET

Electrolyte Analyser



WHOLE BLOOD OR PLASMA
ELECTROLYTE ANALYSER.
ZERO MAINTENANCE.



"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.

i-Smart 30 VET

Electrolyte Analyser

The i-Smart 30 Vet is a whole blood or plasma electrolyte analyser requiring zero maintenance.

FAST & ACCURATE

- Fast, accurate analysis
- Tests electrolytes **Na⁺**, **K⁺** and **Cl⁻**
- 60µl sample volume
- Results in 35 seconds
- Species specific reference ranges
- Direct ISE measurement

ALL-IN-ONE CARTRIDGE

- Cartridge contains electrodes, consumables, waste and sample probe
- No user maintenance
- Automatic calibration

FEATURES & ADVANTAGES

- Large LCD touch screen
- Veterinary specific software
- Battery or mains powered

MEASURED ELECTROLYTES

Analytes	Measuring Range
Na ⁺	20~250 mmol/L
K ⁺	0.5~20.0 mmol/L
Cl ⁻	20~250 mmol/L

PACK SIZES

30, 50, 100, 200 tests

VOLTAGE

100~240 VAC



**Fast, accurate and maintenance free
In-Practice electrolyte analysis**



MEMORY

- 300,000 results storage
- Data list, search and view
- Data downloads to USB
- Print and send capability
- Previous result comparison



SPECIES

- Wide range of species
- Species specific reference ranges reported
- Correlation factor set for each species

+44 (0) 1204 669033 (Option 1)

sales@woodleyequipment.com

www.woodleyequipment.com

Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

**WOODLEY
EQUIPMENT COMPANY LTD.**