

IMMUNO-CHEK G

Foal IgG Test Kit

Immuno-Chek G is an easy to use lateral flow test for semi-quantitative detection of immunoglobulin G (IgG) in a drop of equine whole blood.

Foals are born with no immunoglobulins (IgG) and must receive adequate antibodies from the mare's colostrum (passive transfer of IgG). However, up to 20% of newborn foals experience a partial or total failure of passive transfer of IgG. This can be caused by premature lactation, deficient suckling, malabsorption or low IgG levels in the mare's colostrum. These animals are at high risk of infection and serious illness or death, so it is vitally important that vets undertake routine measurement of the IgG levels of all foals at around 12 hours old to determine there isn't failure of transfer of passive immunity (FTPI).

Quickly identifies immunodeficient foals



Kit Storage

- Room Temperature

Sample

- 10µl EDTA/ Heparin whole blood

Test Time

- 10 minutes


Kit Size

- 5 or 10 Tests

Kit Components

- Pouches containing test device
- Sample diluent
- Buffer
- Pipette & pipette tips

Result Interpretation



Tel: +44 (0)1204 669033

Fax: +44(0)1204 669034

Email: sales@woodleyequipment.com

Web: www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK