

Infusomat Space Infusion Pump

Available for hire or purchase for your clinical trial.



Overview

This Infusomat® Space infusion pump is designed to be smaller, lighter and easier to handle.

This system meets all the demands of an outstanding safety concept. The Infusomat Space is setting new standards in infusion technology.

Infusion rate range - 0.1 - 1200 ml/h.

Operating Features

- Automatic rate calculation by calculating volume over time
- Dose rate mode
- Dose calculation
- Automatic calculation of delivery rate
- Drug Library: up to 1500 drug names can be stored
- Drug specific soft and hard limits
- KOR mode
- Standby mode, adjustable from 1 minute to 24 hours
- Safety Concept:
 - help prompts
 - freeflow protection
 - occlusion alarm
 - bolus reduction
 - upstream pressure detection
 - drop control
 - data lock on 3 safety levels
 - LED alarm indicator with clear message in display
 - Pre - alarms

Specification

- Infusion Rate: 0.1 - 1200 ml/h
- Flow Rate Range: 0.01-1200 ml/h
- Accuracy: \pm 5% over 96 h
- Bolus Delivery Rate: 1 - 1200 ml/h
- Volume Pre-Selection: 0.1 - 9999 ml
- Time Pre-Selection: 1 min - 99 h 59 min
- Occlusion Alarm: 9 occlusion pressure levels (selectable)
- Alarms: Upstream alarm, drop alarm, air alarm, KOR alarm, reminder alarm, battery empty alarm, infusion line not or incorrectly inserted alarm, technical alarm, etc.
- Operating Conditions: + 10°C to + 40°C
- Mains Power Supply: SpaceStation or AC adaptor
- Battery Type: rechargeable NiMH
- Weight: 1.4 kg
- Dimensions: 68h x 214w x 124d mm



Please contact us for a quotation for the hire or purchase of this product

Product description: Infusomat Space Infusion Pump

Product Code: WD7INF011

Version: V3 19/10

Product specifications can be subject to change with



UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company