



# Arterial Blood Gas Syringes

Available in 1ml and 3ml sizes, with or without needle



**SIMPLE, EFFICIENT**  
**BLOOD SAMPLE COLLECTION.**

**"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"**

**WOODLEY**  
EQUIPMENT COMPANY LTD.



# Arterial Blood Gas Syringes

Available in 1ml and 3ml sizes, with or without needle

## FEATURES

- Sterile, single-use device for the collection of whole blood samples
- Luer cap to seal the sample
- Luer lock connection for secure needle attachment
- Breathable plug - double layer porous stone structure, auto-exhaust
- Anticoagulant - ultra-precision sprayed solid lithium heparin, does not dilute blood samples - balanced by calcium ions to ensure accuracy of test results
- Preset system for simple and efficient automatic sampling
- Sharps protection guard on versions with needle to help prevent needlestick injury



## SPECIFICATIONS



Validated for use with  
epoc Critical Care/  
Blood Gas Analyser



### Product Code Description

WD0234	1ml without needle
WD0235	3ml without needle
WD0236	1ml with needle
WD0237	3ml with needle

Supplied in Packs of 50



+44 (0) 1204 669033 (Option 1)



[sales@woodleyequipment.com](mailto:sales@woodleyequipment.com)



[www.woodleyequipment.com](http://www.woodleyequipment.com)



Old Station Park Buildings  
St. John Street  
Horwich  
Bolton  
Lancashire  
BL6 7NY  
UK

**"PASSIONATE ABOUT  
LABORATORY DIAGNOSTICS"**

**WOODLEY**  
EQUIPMENT COMPANY LTD.