

InSight™ DB

Veterinary Chemistry Analyser



- Low cost, All Inclusive Reagent, Support & Maintenance Package
- Includes InSight Lab Manager
- Integrates to your Practice Management System
- Regular Preventative Maintenance Visits by our Nationwide Team of Field Service Technicians
- Out of Hours Technical Support
- Fixed Price for 5 Years



"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.

InSight™ DB Chemistry Analyser

The InSight DB Veterinary Chemistry Analyser is fully automated, uses liquid reagents, and is approximately 70% cheaper to run than dry chemistry.

FEATURES

- User friendly intuitive veterinary software
- Onboard refrigerated reagent/sample compartment
- Minimal user maintenance
- Fast throughput - example batch times:
9 mins 1x GHP, 13 mins 2x GHP, 16 mins 3x GHP
- Continuous loading, no need to wait for batch completion
- 40 reagent positions
- Stat function for urgent tests
- Fully automated with auto-reflex dilution of over range results
- Flexible testing, single parameter or up to 20 user defined profiles
- Online remote assistance directly to your analyser
- Long life ceramic syringes to ensure accuracy and precision
- Environmentally friendly semi-disposable reaction -cuvettes with high efficiency wash system



CUSTOMER DEFINED PROFILES AND/OR TESTS CAN BE PROGRAMMED FROM THE FOLLOWING LIST:

Albumin, Alkaline Phosphatase, ALT, AST, Amylase, Bilirubin (Total), Calcium, Cholesterol, CK, Creatinine, Glucose, GGT, Phosphate, Total Protein and Urea.

BHB, Bile Acids, Chloride, Copper, Fibrinogen, Fructosamine, GLDH, GSH-Px (as Selenium), Haemoglobin, IgG, LDH, Lipase, Magnesium, NEFA, SAA, T4 (Total Thyroxin), Triglycerides, Uric Acid, Urine Creatinine, Zinc, Phenobarbitone, Iron, Ferritin and TIBC

DIMENSIONS/WEIGHT

Dimensions: (W) 414 x (D) 565 x (H) 385 mm

Weight: 28kg

Water Consumption: <3L/H

Available as a low cost package including analyser, reagents and service. Suitable for any size practice, equine practices or for multi-site placements.



Fully supported nationwide by our specialist team of field service technicians.

+44 (0) 1204 669033 (Option 1)

sales@woodleyequipment.com

www.woodleyequipment.com

Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

**WOODLEY
EQUIPMENT COMPANY LTD.**