

InSight[®] DB

Veterinary Chemistry Analyser



**MORE FLEXIBLE & LOWER COST
TO RUN THAN DRY CHEMISTRY**



"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

**WOODLEY
EQUIPMENT COMPANY LTD.**

InSight[®] DB Chemistry Analyser

The InSight DB Veterinary Chemistry Analyser is fully automated, uses liquid reagents, and is **more flexible and cheaper to run than dry chemistry.**

FEATURES

- User friendly intuitive veterinary software
- Onboard refrigerated reagent/sample compartment
- Minimal user maintenance
- Fast throughput - example batch times:
9 mins 1x GHP, 13 mins 2x GHP, 16 mins 3x GHP
- Continuous loading, no need to wait for batch completion
- Stat function for urgent tests
- Fully automated with auto-reflex dilution of over range results
- Flexible testing, single parameter or up to 20 user defined profiles
- Small footprint - 414 (W) x 565 (D) mm
- Online remote assistance directly to your analyser
- Long life ceramic syringes to ensure accuracy and precision
- Environmentally friendly semi-disposable reaction cuvettes with high efficiency wash system



CUSTOMER DEFINED PROFILES AND/OR TESTS CAN BE PROGRAMMED FROM THE FOLLOWING LIST:

Albumin, Alkaline Phosphatase, ALT, AST, Amylase, Bilirubin (Total), Calcium, Cholesterol, CK, Creatinine, CRP, Glucose, GGT, Phosphate, Protein (Total) and Urea.

BHB, Bile Acids, Chloride, Copper, Fibrinogen, Fructosamine, GLDH, GSH-Px (as Selenium), Haemoglobin, IgG, LDH, Lipase, Magnesium, NEFA, SAA, T4 (Total Thyroxin), Triglycerides, Uric Acid, Urine Creatinine & Urine Protein, Zinc, Phenobarbitone, Iron, Ferritin and TIBC.

DIMENSIONS/WEIGHT

Dimensions: 414 (W) x 565 (D) x 385 (H) mm
Weight: 28kg
Water Consumption: <3L/H

**Suitable for any size practice,
equine practices or for
multi-site placements**



+44 (0) 1204 669033 (Option 1)

sales@woodleyequipment.com

www.woodleyequipment.com

Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

**WOODLEY
EQUIPMENT COMPANY LTD.**