

WOODLEY



InSight™ EQA

External Quality Assessment (EQA)

Does your practice use biochemistry or haematology analysers?

External Quality Assessment (EQA) of laboratory equipment is required by the RCVS Practice Standards Scheme and is recommended for any practice with in-practice analysers.

In conjunction with the SRUC (Scotland's Rural College), Woodley Equipment are proud to introduce a 'true' External Assessment with:

InSight™ EQA

A certified independent External Quality Assessment scheme enabling you to provide the best possible patient care.



Why Use EQA?

- Maintaining and improving analytical quality
- Improving inter-laboratory agreement thus raising standards
- Detecting equipment failures and identifying reagent problems
- Reviewing staff training
- Initiating and evaluating corrective actions
- Comparing different analytical methods
- Compliant with the RCVS Practice Standards Scheme

InSight EQA is independent of all analyser manufacturers

To enrol complete the form at the back of this brochure and fax back to 01204 669034, email to sales@woodleyequipment.com or phone Woodley Equipment on 01204 669033 (Option 1)

Tel: 01204 669033 (Option 1)

Fax: 01204 669034

Email: sales@woodleyequipment.com

Web: www.woodleyequipment.com

**WOODLEY
EQUIPMENT COMPANY LTD.**

OLD STATION PARK BUILDINGS, ST. JOHN STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK

Why Choose

InSight™ EQA



- InSight EQA is an internet based clinical chemistry and haematology EQA scheme
- InSight EQA provides a 'true' external assessment by comparing results from multiple analysers of the same make and NOT just with a reference laboratory analyser
- InSight EQA is independent of all analyser manufacturers
- Flexible scheme so you choose which analytes to assess

InSight™ EQA is a **'NOT FOR PROFIT'** scheme in conjunction with **Scotland's Rural College (SRUC)**

SRUC carries out academic and government funded not-for-profit business and is committed to excellence in the advancement, communication and translation of knowledge throughout the rural sector.

It leads on innovation and sustainable development in agriculture, land and the rural sector.

Profits from SRUC's commercial businesses, including External Quality Assessment are used by SRUC to support the not-for-profit businesses, such as SRUC's farming for a better climate which provides practical support to farms to reduce their impact on the climate and the animal behaviour and welfare team who work on understanding behavioural development and create methods of welfare assessment.

For further information on SRUC's activities visit:

www.sruc.ac.uk



To enrol complete the form at the back of this brochure and fax back to 01204 669034, email to sales@woodleyequipment.com or phone Woodley Equipment on 01204 669033 (Option 1)

Tel: 01204 669033 (Option 1)

Fax: 01204 669034

Email: sales@woodleyequipment.com

Web: www.woodleyequipment.com

**WOODLEY
EQUIPMENT COMPANY LTD.**

OLD STATION PARK BUILDINGS, ST. JOHN STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK



Three simple steps

1. We send you a sample to analyse each month
2. You email back the results
3. We send you a report within 7 days comparing your results with all subscribers for both identical analysers and all registered with the scheme

Tailored to suit your practice

InSight EQA offers monthly EQA programmes for **biochemistry**, or **haematology** or both.

You choose what you want to test

You choose which analytes you want to include in your assessment to suit your requirements.

Whatever the species

InSight EQA can cover analytes specific to companion animal or farm species and samples are designed to mimic the performance of a patient sample.

Helping you with quality control

InSight EQA provides an indication of your analyser's accuracy and precision complementing your internal quality control.

What's in the report we send back to you?

A quick reference summary page providing an overview of your performance for each parameter.

- *Number of participants*
- *Mean*
- *Standard deviation (SD)*
- *Standard deviation index (SDI)*
- *Running mean of the SDI*
- *% deviation from the mean*
- *Running mean of the % deviation from the mean*
- *Coefficient of variation*
- *Levey-Jennings and % deviation from the mean charts* allowing quick and easy identification of any trends or bias
- *Histogram* showing the spread of all results

Our reports are sent promptly to your QA Manager so that you can address any problems should they arise.

Certification

We send out certificates at the end of each yearly cycle if your return rate is over 50%.

To enrol complete the form at the back of this brochure and fax back to 01204 669034, email to sales@woodleyequipment.com or phone Woodley Equipment on 01204 669033 (Option 1)

Tel: 01204 669033 (Option 1)

Fax: 01204 669034

Email: sales@woodleyequipment.com

Web: www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS, ST. JOHN STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK



Faxback to:
01204 669034



2017 Prices:

The current standard prices to register for InSight EQA are:

1st April 2017 to 31st March 2018 cycle:

- Single biochemistry analyser: £160+VAT/year
- Single haematology analyser: £160+VAT/year
- One biochemistry and one haematology analyser at the same address: £240+VAT/year
- Analysers at branch surgeries: £100+VAT/year

Prices are reduced pro-rata for practices joining midway through the cycle.

Note that if a practice has a biochemistry analyser and a separate electrolyte analyser at the same address, for charging purposes they are considered as one analyser, although separate reports will be issued for each analyser.

However, if samples are posted to a practice and they choose not to return results that month no credit will be provided.



**Enrol for InSight EQA and receive
a gPet PLUS Veterinary Glucometer Kit,
together with 25 Test Strips
FREE OF CHARGE**

Enrolment Form:

Practice Name: _____ Contact Name: _____

Practice Address: _____ Position/Job Title: _____

_____ Email Address: _____

_____ Tel No: _____

Postcode: _____

On receipt of this request a more detailed Registration Form will be emailed to gather analyser information.

Tel: 01204 669033 (Option 1)
Fax: 01204 669034
Email: sales@woodleyequipment.com
Web: www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS, ST. JOHN STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK