

InSight[®] qLabs[®] PLUS

Veterinary Coagulation Analyser



**COAGULATION SCREEN,
ONE DROP WHOLE BLOOD.
NOW WITH FIBRINOGEN TEST.**

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.

InSight[®] qLabs[®] PLUS

Veterinary Coagulation Analyser

The InSight qLabs Plus Veterinary Coagulation Analyser is calibrated and validated specifically for use in cats and dogs, and allows veterinary practices to carry out point of care coagulation testing in-practice.

FEATURES

- Requires 1 drop of non-anticoagulated whole blood (10µl)
- PT/aPTT tested simultaneously on 1 test strip
- FIB test measures fibrinogen levels within 1-10 minutes from citrated whole blood
- Room temperature strip storage
- Fill test strip directly from syringe
- Liquid QC available
- No need to transfer sample into citrate tube
- 'Test count up' for immediate assessment
- User friendly colour touch screen
- Portable/handheld system with built-in rechargeable battery
- Bluetooth connectivity for software updates
- Charge analyser directly or easily connect to eStation to charge and print
- Automatically print & send results to a PC



User Friendly
Colour Touch Screen

Fibrinogen Test
Now Available



Easily Connect to
eStation to
Charge and Print

SIMPLE TEST PROCEDURE

1. Insert test strip
2. Apply a drop of blood
3. Results displayed on LCD screen

HELPING DIAGNOSE CONDITIONS IN RELATION TO:

- Pre-operative checks for at risk animals
- Toxins and systemic diseases
- Hereditary bleeding disorders
- Excessive bleeding
- Disseminated Intravascular Coagulopathy (DIC)

CARTRIDGE TYPES

- PT/aPTT combo cartridge (12)
- FIB Test Strips (12)

KIT CONTENTS

- qLabs[®] Plus Coagulation Analyser
- Power Adapter
- User Manual
- Lancets (10)
- Carry Case
- Lancing Device
- eStation (optional)

SPECIFICATIONS

Display: 3.5" 320 x 480 colour LCD
Barcode Scanner: QR code scanner
Languages: English, German, French, Italian, Spanish, Portuguese, Polish
Patient Records: Stores up to 2,000 patient records and 500 raw data
Measurable Range:
PT Range: 5-75 seconds
aPTT Range: 6.5-500 seconds
FIB Range: 60-600 seconds

Dimensions: 148 (L) x 70 (W) x 26 (H) mm
Weight: 203g
Operating Conditions:
Temperature: 10°-35°C
Relative Humidity: 10%-90%
Power Supply: 100V-240V AC
Built-in 2500 mAh lithium polymer battery



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

**WOODLEY
EQUIPMENT COMPANY LTD.**