


InSight[®] Roller Mixer

For mixing blood samples, viscous substances
and for liquid-solid suspensions.



GENTLE ROCKING & ROLLING
MOTION FOR OPTIMUM MIXING.




"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.

InSight[®] Roller Mixer

For mixing blood samples, viscous substances
and for liquid-solid suspensions.



FEATURES

- Gentle rocking and rolling motion for optimum mixing
- Fixed speed of 30 RPM with speed accuracy of +/- 2 RPM
- Continuous/infinite mode of operation
- Tilt angle of +/- 3°
- Tool free removal/installation of rollers for replacement and cleaning
- Adjustable rollers spacing to hold different diameter tubes
- Suitable for mixing of blood samples, viscous substances and for liquid-solid suspensions
- Small footprint saving space
- Works with 80% humidity and 40°C
- Protection class of IP 21
- Robust construction for stable operation

SPECIFICATIONS

Motor Type: AC Synchronous Motor
Number of Tube Rollers: 4 Rollers
Rocking Motion: Rocking & Rolling
Speed Range: Fixed 30 RPM
Speed Accuracy: ± 2 RPM
Timer Range: Continuous mode
Tilt Angle: +/- 3°
Max. Loading Capacity: 1.5kg
Removable Rollers: Yes
Roller Size: 330mm (length) x 30mm(roller diameter)
Dimensions: 475 x 210 x 80 mm
Protection Class: IP 21
Permissible Ambient Temperature: 5°C - 40°C
Permissible Relative Moisture: <80%
Weight: 1.5kg
Input Voltage: 110 - 220V 0.02A
Power Consumption: 3W (110V), 6W (220V)



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

**WOODLEY
EQUIPMENT COMPANY LTD.**