

InSight™ V5

Veterinary Haematology Analyser

5 PART DIFFERENTIAL WITH LASER TECHNOLOGY



"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.

FEATURES

- Quick and easy to use
- Small sample volume - 15µl whole blood
- 60 second analysis time
- 5-part differential, 23 parameters
- Species specific reference ranges
- Independent channel for Basophil measurement
- TFT Colour touch screen
- Powerful capability to flag abnormal cells
- Low maintenance
- Fully automated
- Stores up to 50,000 results



SPECIFICATIONS

Principles

- Tri-angle laser scatter method for classifying and counting WBC
- Impedance method for RBC and PLT counting
- Cyanide free reagent and colourimetric method for haemoglobin concentration testing

Parameters

WBC	LYM#	HGB	HCT
NEU%	EOS%	MCV	PLT
NEU#	EOS#	MCH	MPV
MON%	BASO%	MCHC	PDW
MON#	BASO#	RDW-CV	PCT
LYM%	RBC	RDW-SD	

Species :

Dog, Cat, Rabbit, Horse & Cow

User Interface

- Easy to use, menu driven user interface
- Screen Type - LED backlit, colour graphics LCD

Connection Interfaces

- One LAN interface
- Printer Options - Built-in thermal printer or connect to external inkjet and laser printers
- 4 x USB interfaces

Graphic Parameters

- Four 2D scattergrams
- Three histograms

Integrates with Practice Management Systems



Dimensions

- Width 36cm (14.3")
- Depth 43cm (17")
- Height 50cm (19.6")
- Weight - 28kg (61.7lbs)
- Power Supply - 100-240V AC

+44 (0) 1204 669033 (Option 1)

sales@woodleyequipment.com

www.woodleyequipment.com

Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

**WOODLEY
EQUIPMENT COMPANY LTD.**