

InSight V-IA[®]

Veterinary Immunoassay Analyser

An easy to use Immunoassay Analyser providing accurate and reliable results in 3-15 minutes.



**LOW COST,
REFERENCE LAB QUALITY,
IMMUNOASSAY TESTING.**



"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.

An easy to use Immunoassay Analyser providing accurate and reliable results in 3-15 minutes

FEATURES

- Immunofluorescence technology
- Testing time 3-15 minutes (depending on test run)
- Easy to use
- Accurate and reliable
- Touch screen
- Multi-language operation system
- Reagents stored at room temperature
- Maintenance free
- LIS Connection
- Built-in thermal printer



Canine Tests Available:

- cPL
- T4
- TSH
- Cortisol
- Progesterone
- CRP
- D-Dimer
- HbA1C
- NT-proBNP

Feline Tests Available:

- fPL
- T4
- TSH
- Cortisol
- D-Dimer
- SAA
- NT-proBNP

Tests In Development

- Troponin I
- AFP (Alpha-fetoprotein)

SPECIFICATIONS

Dimensions: 213 (W) x 303 (D) x 210 (H) mm
Weight: 4kg
Power Adapter: 100~ 240V, 50~ 60Hz



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

**WOODLEY
EQUIPMENT COMPANY LTD.**