

InSight[®] Vet CGM

Continuous Glucose Monitoring System

CGM now designed and manufactured solely for Veterinarians. Veterinarians and Nurses can now monitor all their Diabetic clients through one easy to use App and Web Portal.



**UP TO 14 DAYS OF
CONTINUOUS MONITORING
FOR DOGS AND CATS**

InSight[®] Vet CGM

Continuous Glucose Monitoring System For Dogs and Cats

Veterinarians and Nurses can now monitor all their Diabetic clients through one easy to use App & web portal

FEATURES

- Comfortable to use
- Accurate results
- Up to 14 days of continuous monitoring
- Calibrated specifically for Dogs and Cats
- Significantly reduces the need for blood sampling, reducing pet's pain, stress and difficulty of monitoring
- Reduced-noise design reduces stress for the pet
- Specially designed software with the unique needs of pets in mind
- App considers factors such as species, size and metabolic rates that influence the pet's glucose levels



CGM now designed & manufactured solely for Veterinarians

INSIGHT VET CGM VS. HUMAN CGM

InSight Vet CGM

- Calibrated specifically for cats and dogs
- Specially designed probe for pets
- The glucose thresholds, parameters and alert ranges are suitable for cats and dogs
- App allows veterinarians account to monitor data for multiple diabetic pets

Human CGM

- Glucose range is designed for humans
- Algorithm is designed for humans
- No after-sales service
- Software is designed for humans

SPECIFICATIONS

Species: Dogs and Cats	Shelf Life: 18 months
Data Interval: 3 minutes	Age: 8 months+
Monitoring: Up to 14 days	Warm-up Time: 30 minutes
Calibration: Factory calibration	MARD: 7.38%
Waterproof: IP27	Prediction Trend: Yes

INSIGHT VET CGM INCLUDES:

Applicator

- Disposable
- Passive activation

Sensor

- Bluetooth connection
- One data point every 3 minutes, monitoring for 14 days
- Flexible probe, comfortable and imperceptible to wear
- IP27 protection rating, water-resistant
- For daily use

App & Web Portal

- Receive data through a compatible smartphone app
- Cloud-based portal for analysing and generating blood glucose reports
- Facilitates sharing with Veterinarians and loved ones



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

WOODLEY
EQUIPMENT COMPANY LTD.