

Jumbo Digital Vaccine Thermometer

Available for hire or purchase for your clinical trial.



Overview

This jumbo digital vaccine thermometer is ideal for use in the clinical trials industry, comes with a triple display and is waterproof.

This thermometer has large digits to aid readability and has several mounting options.

There are 3 different probe types available - please see opposite for further details.

Operating Features

- Triple display: maximum, minimum, current probe temperatures plus room temperature
- 38mm high digits on display
- Sensor and 3m cable perform even when under water
- Alarm alert when temperature rises or falls below a set point, programmable in 1°C increments
- High-impact, chemical resistant ABS plastic case
- Flip-open stand for lab bench, slot for wall mounting, and Velcro™ and magnetic strips to attach to any surface
- Min/max monitors high/low readings
- Supplied with AAA battery and Traceable® Certificate ISO 17025 A2LA to NIST

Specification

- Range: -50°C to +70°C
- Resolution: 0.1°C
- Accuracy: ± 1°C
- Weight: 113g
- Dimensions: 111h x 22w x 8d mm

3 different probe types available

- 1 - Probe sealed in miniature bottle (25 x 64mm) filled with patented nontoxic glycol, bottle insulates sensor from rapid temperature changes when refrigerator door is opened
- 2 - Naked sensor probe measure 4.8mm in diameter, 19mm long
- 3 - 5mL vaccine bottle mimics vaccine temperature, and fits in vaccine tray



Please contact us for a quotation for the hire or purchase of this product

Product description: - Jumbo Digital Vaccine Thermometer

Product Code: -

Product specifications can be subject to change without prior notice. Version: V2 10/19

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company