



Veterinary Lactate Meter

Easy to use handheld meter providing
blood lactate results in just 45 seconds



FAST, ACCURATE
LACTATE RESULTS.
5 YEAR WARRANTY.

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.

L-Pet[®] Veterinary Lactate Meter

An easy to use handheld Lactate Meter providing blood lactate results in just 45 seconds.

L-Pet only requires a small drop of blood (<3µl).

The test strips draw blood in and then measure using an advanced biosensor specific for lactate concentration in blood.

KEY FEATURES

- Advanced Biosensor Test Strips
- Small sample size (<3µl)
- Results in less than 45 seconds
- Large LCD display
- Stores up to 300 results
- Easy operation
- Portable
- Specific Lactate Analyser QC Material available



Easy to use handheld meter providing blood lactate results in just 45 seconds

SPECIFICATIONS

Measuring Range: 0.7 ~ 22.2 mmol/L (6-200 mg/dL)

Sample Volume: 3µl whole blood

Measuring Time: 45 seconds

Operating Environment: 0°C ~ 40°C (50°F ~ 104°F)
Humidity 20% ~ 80% (RH)

Battery Life: Approximately 1000 tests

Calibration: External Code Card

Memory Capacity: 300 sets with date and time

Temperature Compensation: Auto-sensing compensation

Battery Source: 1x CR2032 3.0V Lithium Battery

PC Interface: RS 232

Audible Beep: Yes

Weight: 64g with battery

KIT CONTENTS

- L-Pet Meter
- L-Pet Test Strips (5)
- Lancing Device
- Lancets (5)
- Monitor Checker
- Battery
- Code Card
- Wallet
- User Guide

+44 (0) 1204 669033 (Option 1)

sales@woodleyequipment.com

www.woodleyequipment.com

Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

**WOODLEY
EQUIPMENT COMPANY LTD.**