

Labofuge® 200 Centrifuge 110V

Available for hire or purchase for your clinical trial.



Overview

The functional Labofuge® 200 is ideal for use in medical practices, clinical and small laboratories and as a stand-by unit in large laboratories. Offered with rotor and all necessary adapters, the Labofuge® 200 centrifuge can be run the day it is received.

The Labofuge® 200 is equipped with a brushless motor, eliminating the need to replace carbon brushes. In addition to time and cost savings, this also ensures clean operating conditions.

Operating Feature

- Easy to use
- Bright digital displays and touch-pad keys offer quick and easy setting of speed and run-time
- START key repeats runs
- User-friendly self-diagnostics indicate errors on the display
- Brushless motor is virtually maintenance-free and ensures clean operating conditions
- Made of impact resistant, fiberglass reinforced polyamide
- Equipped with lid lock, lid interlock and steel armored guard ring
- Bright digital displays and touch-pad keys
- Supplied adapters permit centrifugation of all standard 5, 7, 10 and 15mL tubes, including the popular Monovette, Vacutainer and VenoJect blood collection tubes
- Includes: 12 x 15mL fixed-angle rotor and adapter set; adapter set includes 12 each of 7mL (yellow), 5mL (cream), and rubber distance pads

Specification

- Max. Speed: 5300 RPM
- Min. Speed: 1600 RPM
- Max. RCF x g: 3030
- Min. RCF x g: 270
- Max. Capacity ml: 12 x 15 (glass)
- Controller: microprocessor controller for time and speed
- Runtime: 1 - 99 min. and continuous operation
- Program Memory: stores values last entered, unlimited
- Power Consumption: 65 W
- Weight: 10.7 kg
- Dimensions: 240h x 284w x 375d mm



Please contact us for a quotation for the hire or purchase of this product

Product description: Labofuge® 200 Centrifuge

Product Code: WD7CEN045

Product specifications can be subject to change without prior notice

Version: V2 10/19.

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company