

Laboratory Undercounter Freezer, 151 L



Available for hire or purchase for your clinical trial.

Overview

This solid door 151 litre undercounter laboratory freezer is spark free and ideal for the safe storage of a wide variety of liquids and chemicals.

With new unique Anti Microbial Handle to help prevent the growth of bacteria and minimise cross contamination.

The external digital controller records temperatures to an accuracy of 0.1°C and provides precise temperature readings without having to open the door. This prevents unnecessary temperature fluctuations within the freezer.

Operating Features

- Solid, hinged, reversible door
- Open door alarm
- Lock
- 2 adjustable feet, 2 rear rollers
- Temperature display in the panel
- Temperature alarm
- 3 wire drawers
- Manual freezer defrost
- Colour: White/Green

Specification

- Gross Capacity: 151 Litre
- Net capacity: 134 litres
- Operating Temperature: - 18°C to - 25°C
- Internal Finish: White
- Door Type: Hinged, reversible
- Door Style: Solid
- Spark Free
- Climate Class: SN
- Voltage Supply: 230 V / 50 Hz (Only)
- Weight: 36 kg
- Exterior Dimensions: 845 x 595 x 595 mm HxWxD



Please contact us for a quotation for the hire or purchase of this product

Product description: Undercounter Laboratory Freezer 151 Litre

Product code - WD7FRZ006

Version: V2 10/19

Product specifications can be subject to change without prior notice.

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company