

# Laboratory Freezer, 50 Litre

Available for hire or purchase for your clinical trial.



## Overview

The -25°C benchtop laboratory freezer has an integral control panel on the outside of the fridge which sits flush with the door and its easy to use design makes checking the internal temperature straightforward.

As there is no need to open the door to check the temperature it keeps the internal air temperature more stable. The spark-free interior reduces the risk of internal explosion.

Come complete with door locks and two keys ensuring that only authorised individuals are able to access them.

## Operating Features

- Solid, hinged, reversible door
- Lock
- Fixed feet
- Temperature display in the panel
- Temperature alarm
- 1 wire shelf
- Manual freezer defrost

## Specification

- Colour: white
- Net capacity: 38 litres
- Gross capacity: 50 litres
- Door Type: hinged, reversible
- Door Style: solid
- Operating Temperature: - 18°C to - 25°C
- Climate Class: N
- Voltage Supply: 230 V / 50 Hz
- Weight: 17 kg
- Exterior Dimensions: 562h x 500w x 510d mm



**Please contact us for a quotation for the hire or purchase of this product**

Product description: - Laboratory Freezer, 50 Litre

Product Code: WD7FRZ004

Version: V3 10/19

Product specifications can be subject to change without prior notice.

### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: [hello@woodleytrialsolutions.com](mailto:hello@woodleytrialsolutions.com)

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company