

Laboratory Refrigerator, 107 Litre

Available for hire or purchase for your clinical trial.



Overview

This 107 litre under counter medical/laboratory refrigerator has an integral control panel on the outside which sits flush with the door. Its easy to use design makes checking the internal temperature straightforward. As there is no need to open the door to check the temperature, it keeps the internal air temperature more stable. It also features automatic defrost to ensure more efficient cooling.

The spark-free interior reduces the risk of internal explosion. Security is essential in laboratories so this fridge come complete with door locks and two keys.

Operating Feature

- 2°C to 10°C temperature range
- Spark free interior
- External digital controller and temperature displays
- Fitted lock and 2 keys
- Audible and visual temperature alarms

Specification

- Colour: White
- Gross Capacity: 107 Litres
- Nett Capacity: 103 Litres
- Temperature range: 2°C to 10°C
- Door Type: Solid
- Door Style: Hinged reversible
- Climate Class: N
- Refrigeration Type: Forced air refrigeration system
- Voltage Supply: 230 V
-
- External Size: 495 x 620 x 830 mm (WDH)
- Shipping Dimensions: 560 x 850 x 655 mm (WDH)
- Gross Weight: 31.5kg



Please contact us for a quotation for the hire or purchase of this product

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: enquiries@woodleytrialsolutions.com

www.woodleytrialsolutions.com

Product description: Laboratory Refrigerator

Product Code: WD7FRG373

Product specifications can be subject to change without prior notice.

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a divisions of Woodley Equipment Company