

# Laboratory Refrigerator 151 Litre

Available for hire or purchase for your clinical trial.



## Overview

This 151 litre under counter medical laboratory refrigerator is ideal for the safe storage of a wide variety of volatile liquids and chemicals. Advanced controller records temperatures accurate to 0.1° C.

The unit has audible and visual high and low temperature alarms.

The spark-free interior reduces the risk of internal explosion. Security is essential in laboratories so this fridge comes complete with door locks and two keys.

## Operating Features

- Operates at +6° C for best possible storage of samples and chemicals
- Lock and two keys
- External digital controller for temperature readings without opening the door so preventing unnecessary temperature fluctuations within the fridge.
- Solid, reversible, hinged door
- Anti microbial handle
- 3 wire shelves with six positions
- 2 wire adjustable door bins
- Automatic fridge defrost

## Specification

- Colour: White / Green
- Net capacity: 145 litres
- Gross capacity: 151 litres
- Door Type: hinged, reversible
- Door Style: solid
- Operating Temperature: + 2°C to +10°C
- Climate Class: SN
- Voltage Supply: 230 V / 50 Hz
- Weight: 36kg
- Exterior Dimensions: 845 x 595 x 595mm H x W x D



**Please contact us for a quotation for the hire or purchase of this product**

Product description: - Laboratory Refrigerator, 151 Litre

Product Code: WD7FRG018

Version: V2 10/19

Product specifications can be subject to change without prior notice.

### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: [hello@woodleytrialsolutions.com](mailto:hello@woodleytrialsolutions.com)

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company