

Laboratory Refrigerator, 82 Litre

Available for hire or purchase for your clinical trial.



Overview

This 82 litre benchtop medical/laboratory refrigerator has an integral control panel on the outside which sits flush with the door. Its easy to use design makes checking the internal temperature straightforward. As there is no need to open the door to check the temperature, it keeps the internal air temperature more stable. It also features automatic defrost to ensure more efficient cooling.

The spark-free interior reduces the risk of internal explosion. Security is essential in laboratories so this fridge come complete with door locks and two keys.

Operating Features

- Open door sensor
- Lock
- Fixed feet
- Temperature display in the panel
- Temperature alarm
- 2 shelves
- 1 wire drawer
- Interior light
- Automatic fridge defrost

Specification

- Colour: white
- Net capacity: 67 litres
- Gross capacity: 82 litres
- Door Type: hinged, reversible
- Door Style: solid
- Operating Temperature: + 2°C to +10°C
- Climate Class: N
- Voltage Supply: 230 V / 50 Hz
- Weight: 27.5 kg
- Exterior Dimensions: 660 x 502 x 540 mm HxWxD



Please contact us for a quotation for the hire or purchase of this product

Product description: Laboratory Refrigerator, 82 Litre

Product Code: WD7FRG016

Version: V2 10/19

Product specifications can be subject to change without prior notice.

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company