

# StatStrip® Lactate Xpress Meter

Results in 13 seconds



**FAST, ACCURATE**  
**LACTATE RESULTS.**

**IN 13 SECONDS.**

*nova*  
biomedical



**"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"**

**WOODLEY**  
EQUIPMENT COMPANY LTD.

# StatStrip® Lactate Xpress Meter

StatStrip Lactate Xpress Meter for use in the quantitative determination of Lactate in venous and arterial whole blood.

A hand-held, battery-powered meter to measure lactate in a whole blood sample.

The meter can store up to 400 patient and quality control test results.

## KEY FEATURES

- Measures lactate in whole blood
- Easy to use
- Portable
- No calibration required
- Accurate, quantitative results in 13 seconds
- Small 0.6µl sample size
- Venous or arterial whole blood
- Stores 400 results
- Corrects results for HCT
- Eliminates interfering substances



## LACTATE MEASUREMENT RANGE

0.3-2.0 mmol/L  
2.7-180 mg/dL

## SAMPLE TYPES

Whole Blood: Venous & Arterial

## SPECIFICATIONS

Temperature: 15-40°C (59-104°F)  
Weight: 75g (2.65 oz)  
Size: 91.4 x 58.4 x 22.9 mm

*nova*  
biomedical

+44 (0) 1204 669033 (Option 1)

sales@woodleyequipment.com

www.woodleyequipment.com

Old Station Park Buildings  
St. John Street  
Horwich  
Bolton  
Lancashire  
BL6 7NY  
UK

**"PASSIONATE ABOUT  
LABORATORY DIAGNOSTICS"**

**WOODLEY**  
EQUIPMENT COMPANY LTD.