

# M24/7 Ambulatory BP Monitor

Available for hire or purchase for your clinical trial.



## Overview

The M24/7 Ambulatory Blood Pressure Monitor is fully automatic blood pressure monitor. Windows based software is included for easy set up and analysis of results.

Extensive testing by independent clinicians have validated the accuracy of this ambulatory blood pressure monitor.

The monitor offers doctors and investigators the opportunity to complement their current diagnostic methods by monitoring patients away from the clinic for up to 48 hours

## Operating Features

- Fully automatic blood pressure monitor
- Blood pressure maximum storage: Over 400 automatic measurements
- Blood pressure measurement method: Oscillometric
- Liquid crystal display
- Compact, lightweight design
- Data transmission: USB optical cable
- PC Interface: USB interface
- Accessories available
- 2 cuffs supplied
- Packaged in a soft pouch
- Dimensions 70(h) x 99(w) x 30(d) mm
- Weight approximately 240g including batteries
- Power source: Four 1.5V AA batterie

**Please contact us for a quotation for the hire or purchase of this product**

### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: [hello@woodleytrialsolutions.com](mailto:hello@woodleytrialsolutions.com)

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

## Specification

- Operating Environment
  - Temperature: 10 to 45.0°C
  - Humidity: 10 - 95% non condensing
  - Atmospheric pressure: 70 - 160 KPa
- Pressure Measurement Range: 0 - 300mmHg
- Status Accuracy:  $\pm 3$ mmHg
- Pulse Rate: Measurement Range: 40 - 200 BPM
- Pressure Sensor: Piezo -resistive
- Safety: Maximum inflation 300mmHg safety release valve
- Deflation and Rapid Air Release: Automatic pressure release valve.



Product description: 24/7 Ambulatory Blood Pressure Monitor  
Product Code: WD7BPM016



Product specifications can be subject to change without prior notice

Version: V2 10/19

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company