

MAC 1600 ECG Machine

Available for hire or purchase for your clinical trial.



Overview

The MAC1600 ECG machine can be easily configured to meet the needs of nearly every environment. With instant assessment of the ECG signal quality and an easy connection to your ECG cardiology management system, the new MAC1600 could not be smarter.

The MAC 1600 offers more than just acquisition; it's online editing and re-analysis of ECG waveform capabilities enhance the accuracy of patients' records.

Operating Features

- Soft key menu
- 7" full-colour digital display
- Automatically captures 10-second ECG with one-touch operation
- Rhythm ECG mode to continuously review all 12 leads
- Arrhythmia mode for an extended review of arrhythmias or events
- Integrated handle making transportation easy
- Slide out paper drawer
- Paper Speed: 5, 12.5, 25 or 50 mm/sec
- Online editing and re-analysis
- 6 or 12 lead display
- Automatic storage, transmission and deletion of data
- Barcode reader (optional)
- Save up to 100 ECG records (optional)
- Supports stress testing with treadmill, bicycle or master step stress devices (optional)

Specification

- Display: 7" full colour digital 800 x 480 pixels
- Power Supply: 100 to 240V, +10-15%, 50-60Hz \pm 10%
- Battery Operation: user replaceable 16.4V 2300AH Lithium Ion
- Battery Capacity: 50 single page reports (typical) 3 hrs continuous display (without printing)
- Battery Charge Time: approximately 2.5 hrs from total discharge
- Operating Temperature: 10°C to 40°C
- Storage Temperature: -40°C to 70°C
- Operating Humidity: 20 - 95% RH
- Storage Humidity: 15 - 95% RH
- Operating Pressure: 700 to 1060 hPa
- Storage Pressure: 500 to 1060 hPa
- Weight: 6kgs with battery installed
- Dimensions: 400h x 350w x 150d mm



Please contact us for a quotation for the hire or purchase of this product

Product description: MAC1600 ECG Machine

Product Code: WD7ECG026

Version: V2 10/19

Product specifications can be subject to change without prior notice



UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a divisions of Woodley Equipment Company