

Mechanical Flat Scales

Available for hire or purchase for your clinical trial.



Overview

These mechanical flat scales with 1 lbs graduation provide precision, reliability and incorporate a good design.

They feature a large, clear and well defined dial making it easy to read the measurement.

They have a non-slip platform making it easy for the patient to step on and off.

Operating Features

- Dial:
 - large area
 - well-defined
- Platform:
 - large area
 - low level
 - non-slip
 - scratch resistant
- Weight display immediately after stepping on the scale
- Rugged metal housing withstands longterm use
- Galvanised frame, load lever, counterweight and all sheet metal parts

Specification

- Capacity: 320lbs/150 kg
- Graduation: 1 lbs, 500g
- Functions: reset-to-zero function
- Weight: 7.7 lbs / 3.5kg
- Dimensions (WxHxD): 12 x 4.7 x 18.5" / 305 x 120 x 470 mm



Please contact us for a quotation for the hire or purchase of this product

Product description: - Mechanical Flat Scales

Product code - WD7PAT089

Product specifications can be subject to change without prior notice. Version: V2 10/19

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company